

SERIE 4 - RAPID, RELIABLE, EFFICIENT





www.ecolab.com

www.soluscope.com

The **power** at your **fingertips**

Using cutting edge technologies for improved performances.

WORKFLOW OPTIMISATION

The Soluscope Serie 4 PA offers fast reprocessing cycles to adapt your organization to growing endoscopy activity.

A UNIQUE AND PATENTED DUAL WASH TECHNOLOGY

Semi-immersion bathes the endoscope, whilst the shower adds mechanical cleaning action.

Better results using less water and chemistry for reduced global costs.

LIQUID CHEMICAL STERILANT (LCS) CYCLE

The LCS cycle achieves higher level of efficacy against bacterial endospores (compared to high level disinfection) for challenging cases. Cycle efficacy has been tested based on ISO 14937 protocol.







User-centered **design**

Advanced ergonomics for fast adoption.

SMART AND EASY-TO-USE INTERFACE

The Soluscope Serie 4 PA brings the benefit of an intuitive and graphic interface for fast adoption and quick use. The large and comfortable 7'' touch screen enables easy access to all cycles and machine parameters.

AN INTEGRATED CONCEPT

The Soluscope Serie 4 PA is easy to install. All the elements (chemistry, filters, traceability, barcode reader) are integrated and allow easy access for maintenance.





TURNKEY ELECTRONIC TRACEABILITY AND DATA MANAGEMENT

The Soluscope Serie 4 PA includes barcode reader, printer and local data storage.

It is compatible with our EP-Digital suite for electronic traceability, real-time monitoring of endoscope inventory and enables microbiological samplings and services records.



MINIMISED HANDLING TO REDUCE CROSS-CONTAMINATION RISKS

Dedicated connectors reduce the need for endoscope handling and subsequently the risk of contamination.



A RELIABLE AND REPRODUCIBLE PROCESS

A double check system with alarm verifies all critical parameters to ensure constant performance.

A continuous leak test checks sheath integrity to guarantee endoscope safety during the whole process.

A patented volumetric channel control monitors each individual channel.

Highest standards to ensure **patient safety**

CE MARK / ISO 15883 -1 & 4 compliant.

PROVEN SAFETY

Proven and tested according to EN ISO 15883-4:2018, Soluscope Serie 4 PA in combination with its dedicated, specifically designed chemistry, is compatible with a large range of flexible endoscopes.

CONNECTORS

Compatible with all Soluscope equipments, the connector stays hooked to the endoscope during the entire reprocessing journey, from pre-cleaning to storage.

This link ensures both connector and endoscope are at the same level of hygiene and provides an efficient and safe process before use on a patient.







SOLUSCOPE PAA

- Peracetic acid-based disinfectant
- Broad spectrum antimicrobial activity
- Compatible with a large range of flexible endoscopes

SOLUSCOPE CLN

- pH neutral detergent
- Proven efficacy against biofilm and proteins
- Non foaming

"Ecolab provides a complete solution for flexible endoscope reprocessing"

SERI

SERIE PA



THE ECOLAB SOLUTION

Ecolab provides a complete solution for flexible endoscope reprocessing.

It offers products for every step of the reprocessing cycle.

Each part is certified, easy to use and effective in order to reduce your daily load.



Non-contractual pictures - Soluscope reserves the right to make any changes to the specifications given in this document without notice or obligation.



SOLUSCOPE SAS (An Ecolab Company) 100, rue du Fauge Z.I Les Paluds - 13400 Aubagne - France Distributed by: **ECOLAB EUROPE** Richtistrasse 7 8304 Wallisellen Switzerland





www.soluscope.com

www.ecolab.com