SOLUSCOPE TRANSPORT

the next generation of storage & transport solution



SOLUSCOPE TRANSPORT

Flexible



Multiple-use devices

ENHANCED MANAGEMENT OF YOUR ENDOSCOPES

The Soluscope Trolley, Soluscope Rack and Soluscope EasyTag are designed to optimise the endoscope turnaround between examination and reprocessing rooms. They are compatible with DSC8000 storage cabinet trays and Soluscope transport trays.

SOLUTION FOR ALL SITUATIONS

Beyond its transport function, the *Soluscope Trolley* is suitable to store up to 4 ready-to-use endoscopes for multiple exams or when the examination room is at a distance from the reprocessing room. It can also be used as a cart at the patient's bedside thanks to its ideally positioned tabletop.

SOLUTION FOR ALL ORGANIZATIONS

Soluscope EasyTag labels are made of 2 self-adhesive green and red stickers allowing you to quickly see if your endoscope is ready to be used or if it has to be reprocessed. Intuitive, they are easily and transparently integrated into the daily organization of all endoscopy units.

FOR ALL KINDS OF ENDOSCOPES

The *Soluscope Rack* temporarily stores clean endoscopes or rarely used endoscopes. It can also store dirty endoscopes after exams while awaiting reprocessing.

Ergonomic



be Trolley & Rack offers a concept specially designed for endoscope handling and easy.

This solution improves your team's daily workflow in an efficient manner

giving you comfort and peace of mind.

Your life made easier

MODULARITY

The Soluscope Trolley is compact and adapted to small spaces. It can transport up to 4 endoscopes at any given time. It includes three accessories you can use according to your needs: a tabletop for surgical drapes, sterile gloves or other medical accessories, a document holder to transport patient files, and a special compartment to transport cleaning agents used to disinfect the tray. Its large noiseless wheels make it a pleasure to handle, not to mention its outstanding original and modern design.

OPTIMIZATION

Large in size and easy to remove, *Soluscope EasyTag* improves the visible state of your endoscopes allowing better efficiency in managing your daily endoscope flow.

10 ENDOSCOPES

The *Soluscope Rack* can store up to 10 trays. Its small footprint can allow you to fit it in small spaces. It's the ideal accessory to *Soluscope* AERs and cabinet portfolio to extend your storage capacity.

EASY MAINTENANCE

The Soluscope Trolley and the Soluscope Rack are easy to clean with a surface detergent or disinfectant solutions.

Safe



HIGH-LEVEL PROTECTION

Once positioned in the tray and covered with its individual lid, the endoscope is protected from physical shocks or kinking and from cross-contamination during transportation to the examination room. Medical staff are protected from any microbiological contamination.

Moreover, the Soluscope Trolley and the Soluscope Rack are equipped with sliding rails to prevent the trays from falling.

A CLEAR VISION

Wherever your endoscope is, Soluscope EasyTag indicates if your endoscope is ready to use, avoiding errors and keeping patients safe. The visible status of the endoscope makes it possible to protect the user from all risks of contamination.

STABLE

The guiding wheels of the *Soluscope Trolley* are equipped with brakes providing better stability.



TECHNICAL SPECIFICATIONS - TROLLEY				
Capacity	4 trays			
Dimensions (W x D x H)	50 x 60 x 87 cm (with handles) 50 x 60 x 92 cm (with a tray on top shelf)			
Weight empty	20 kg			


Capacity	10 trays
Dimensions (W x D x H)	61 x 61 x 193 cm
Weight empty	40 kg
	_

TECHNICAL SPECIFICATIONS - TRANSPORT TRAY				
Capacity	1 endoscope			
Dimensions (W x D x H)	55 x 48 x 15 cm			
Weight empty / with lid	1.5 / 2.1 kg			

TECHNICAL SPECIFICATIONS – EASYTAG			
Capacity	500 labels per roll		
Size	10.4 cm x 10.4 cm		







High Level Disinfection

SOLUSCOPE COMPLETE SOLUTION

Soluscope proposes a complete solution for flexible endoscope reprocessing.

It offers products for every step of the reprocessing cycle.

Each part is certified, easy to use and effective in order to reduce your daily load.



www.soluscope.com

Soluscope products are compliant with European and international standards, and CE marked. ISO 9001 version 2015 and ISO 13485 version 2016 certified, Soluscope is continuously committed to satisfying its clients rigorous demands.

a	istr	ΊΟL	itea	bу

31638-FN-I F